

Experience the Healing Nature of Light

1 + 1 = 3

Acupuncture

Low Level Laser Therapy

*Expand Your Practice
Enhance Clinical Outcomes
Increase Revenue & ROI*

The MedX Acupuncture Laser is an ergonomically designed, light weight and comfortable to use. Stimulate meridian points for all of the benefits of acupuncture - without the needles. With the added positive benefit of low level laser therapy, a proven, drug free solution to eliminate pain, decrease inflammation, alleviate swelling and accelerate tissue repair. Additionally, laser acupuncture is an effective treatment for smoking cessation and weight management.

In addition to laser acupuncture, the MedX Portable Laser can treat directly over injured lesions, tissue, nerve roots, trigger points and referred areas of pain. Rapid pain relief, substantial reduction in inflammation, swelling and accelerated tissue repair - accomplished without needles. See the MedX Treatment Manual for over 65 protocols. The treatments incorporate laser stimulation of the injured tissue or lesion, referred areas of pain, nerve roots and acupuncture points and have a prolonged treatment effect.

The merging of acupuncture's ancient wisdom and research **plus** the more recent, high tech proven low level laser combines the best of both worlds.

- ↪ Expand your practice by treating a greater variety of condition.
- ↪ Enhance your clinical results with faster, more reliable patient outcomes.
- ↪ Increase your revenue and quickly pay off your equipment with increased patient referrals, additional reimbursement categories and more time-efficient treatment protocols.

Combine the ancient art of acupuncture and the innovative laser technology. Acupuncture achieves systemic results and rebalancing of qi for improved health. Laser therapy provides more localized, focused outcomes from direct energy stimulation over the affected tissue or source of pain. Together these two approaches produce greater outcomes than either can provide separately.



**MedX Acupuncture
Portable Laser**

See how the MedX Acupuncture Laser can expand your practice, augment your clinical outcomes and boost your reimbursement and revenue.





EXPERIENCE THE BENEFITS OF LOW LEVEL LASER THERAPY... AT THE SPEED OF LIGHT

Now there is a treatment for smoking cessation and weight management. By stimulating the body's production of endorphins and serotonin, the cravings for cigarettes and food are diminished while the behaviour is altered.

Benefits of Low Level Laser Therapy

- ↻ Eliminate or significantly reduce pain
- ↻ Accelerate soft and hard tissue repair
- ↻ Decrease inflammation
- ↻ Minimize swelling
- ↻ Relax muscle tension
- ↻ Resolve muscle spasms

Ask about our Smoking Cessation and Weight Management Programs.....



MedX Acupuncture Portable Laser

Category	Description
Laser Diodes	GaAlAs
Laser Equipment Safety Classification	Class IIIb
FDA Regulatory Classification	Class II
Regulatory Approvals	FDA Cleared, Health Canada Licensed, CE Marked, ISO 13485:2003 Registered
Nominal Laser Wavelength	808 nm (near-infrared)
Guide Light	633 nm (visible red)
Duty Cycle	100% (Continuous mode)
Total Optical Output Power	200 mW +/- 20%
Max. Optical Output at Lens.	250 mW +/- 20% (without light guide)
Output	Continuous: 2 Joules - 10 Seconds 4 Joules - 20 Seconds 8 Joules - 40 Seconds
Laser Power Monitoring	Automatic, continuous
Max Beam Divergence	17 x 40 degrees (prior to lens) (Collimated Output)
Spot size	10 mm ² +/- 20%
Safety Interlock	Yes (Remote Safety Interlock Optional)
NOHD (Nominal Ocular Hazard Distance)	2.6 meters (worst case) 33 cm (typical)
Battery	3.6 V _{DC} Li-Ion, 1500 mAh capacity
Dimensions	-23cm x 8cm x 4cm
Weight	-190g (not including accessories)

The Acupuncture Portable Laser kit includes:

- 1 Acupuncture Portable (LPS 200) Laser System
- 1 Recharger: 9V_{DC}
- 1 Acupuncture Light Guide with metal sleeve
- 1 Safety Key
- 1 Padded Carrying Case
- 2 Standard Protective Goggles MPG01 (one each for patient and clinician)
- 1 Operating Guide
- 1 MedX Treatment Manual (with over 65 musculoskeletal protocols)
- 1 MedX Auricular Acupuncture Points Charts for Smoking and Weight Management

