




## Light / Laser Therapy Products



The **MedX Rehab Console System** has been developed to offer product flexibility and choice. The components have been designed for ease of use by licensed healthcare professionals and/or their delegated assistants.

The **MedX Console** powers the convenient, high-quality laser and light therapy accessories, including both a low level laser therapy device of choice and two super luminous diodes (SLD) accessories.

The **MedX Rehab Laser accessory** is a handheld low level laser therapy device for deeper tissue conditions. This accessory features three near-infrared laser diodes for therapeutic applications and a visible LED guide light.

The **MedX Dental** and **Acupuncture Laser accessories** are also handheld and provide low level laser therapy to a smaller focused treatment area or acupuncture point. This single near infrared laser diode tip has a visible LED guide light and can be used with either the intra-oral or extra-oral dental or acupuncture light guides.

The **MedX SLD Cluster accessory** offers SLD laser and light therapy for broader applications. This accessory can be strapped to the body for "hands-free", unattended treatment.

The **MedX<sup>HOME</sup> System** has been developed to offer a home care option to patients, while being monitored by the clinician. The components have been designed for ease of use by licensed healthcare professionals and their patients.

The lightweight, **MedX<sup>HOME</sup>** unit controls the MedX SLD light therapy cluster. The MedX SLD cluster uses SLD technology in both infrared and visible red diodes (52 diodes) which produce photons of energy and trigger physiological changes in tissue. With the included 18" flexstrap, the MedX SLD cluster can be strapped to the body and be used hands-free. Or, if preferred, it can be comfortably handheld over the area to be treated.



**MedX** has 3 **Portable Lasers (Rehab, Oralase and Acupuncture)** which deliver high-quality laser and light therapy with the convenience of portability. The ergonomic design makes the device functional, easy to use and comfortable to operate.

For the Oralase and Acupuncture Portables MedX has made available different size tips for even more versatility.