

EXPERIENCE THE HEALING NATURE OF LIGHT



MEDX PHOTOTHERAPY 200mW PORTABLE LASER

MedX® Phototherapy is a proven and highly effective treatment which increases local microcirculation: accelerates tissue repair, reduces pain, decreases inflammation/edema, and enhances immune response and angiogenesis.

The MedX low level laser delivers all the power of the corded device with the added flexibility of remote treatment while on the go. Automatic, timed cycles with audible feedback of energy delivered make treatments simple and reliable. The unit is ergonomically designed and is focused for deeper tissue penetration. Three laser diode cluster configuration treats larger areas than a single diode laser. The MedX portable laser has three infra red GaAlAs laser diodes, plus a visible guide LED. The unit is lightweight with a built in rechargeable NiMH nickel metal hydride battery.

The MedX Treatment Manual illustrates comprehensive protocols for a wide variety of common conditions, enabling the clinician to effectively treat patients to accelerate healing and reduce pain for improved outcomes.


Phototherapy has gone through stringent regulatory processes. MedX products are US FDA and Health Canada cleared and are produced in an ISO 13485: 2003 certified manufacturing and testing facility. For a complete profile of MedX Health Corp. and its products visit www.medxhealth.com

MedX phototherapy stimulates the body's natural repair response at the cellular level to reduce pain and accelerate healing.




1.888.363-3112 (US/Cdn only) • 905.670-4428 • info@medxhealth.com • www.medxhealth.com

EXPERIENCE THE HEALING NATURE OF LIGHT



MEDX PORTABLE LASER KIT

INCLUDES:



©MedX Health Corp. 2005

Hand-held low-level laser:

The MedX laser provides a total of 200mW of power with three GaAlAs laser diodes and a visible guide LED. Class IIIb laser.

Carrying case

MedX Treatment Manual

Highly effective and efficient protocols illustrate how to deliver SLD and laser treatments for over 65 conditions.

Operating manual:

Two sets of goggles:

Patient and clinician goggles.

Power supply

9V Battery Charger

MedX Kits: Product Description

0904	MedX 200mW Portable Laser
1035	MedX ^{HOME} System
1116	MedX Console + 200mW Laser + 2 SLDs

